

SCINDERM MEDICAL
"Evidence Based Skin Care"



KojiDerm Brightening Pads

w/SKIN LIGHTENING ADDITIVE
High Potency
Skin Brightening System




Prescribing Information

KojiDerm-CAB Brightening Pads

- Contains enzyme-activated Arbutin for enhanced skin lightening and brightening
- Also contains penetration-enhancing Salicylic acid
- Fortified with Kojic and Ascorbic acids for chelating and anti-oxidant benefits
- Compounded freshly on dispensing for optimal potency
- Lightweight fragrance—and preservative-free solution vehicle

KojiDerm-CAB products contain enzyme activated Arbutin in combination with Kojic, Ascorbic and Salicylic acids. Derived from the Bearberry plant, Arbutin is a chemically classified as a phenolic glucoside having the following chemical formula and structure:



KojiLac-CAB products are intended for all skin types except sensitive

Directions: Apply 1-2 times, daily or as directed.

KojiDerm-CAB Brightening Pads:

For face and body skin lightening and brightening. Use twice daily.

How supplied: KojiDerm-CAB Pads: 60 Pads

With pad solution and vial containing skin lightening additives.

To assure maximum potency, The pads is supplied as a solution vehicle accompanied by a sealed vial containing a mixture of Arbutin, Beta glucosidase, Kojic and Ascorbic acids.

At the time of dispensing, the powder is incorporated into the solution and is stable for 60 days when stored at controlled room temperature. Supplied with a sticker to mark the day of expiration.

Contains (on compounding): Water, Ethoxydiglycol, 10% SD Alcohol 40-B, Ascorbic acid, Pentylene glycol, Arbutin, Kojic acid, Salicylic acid, Butylene glycol, Camellia oleifera extract, L-Ergothioneine, Evodia rutaecarpa fruit extract, Glycerin, Rumex crispus extract, PEG-8 Dimethicone, Laureth-4, Trolamine, Sodium metabisulfate, Edetate disodium and Beta glucosidase.

SKIN TIP>>> For stubborn hyperpigmentation/melasma. Best results combine with Kojiderm Cream with Arbutase for stubborn melasma/chloasma.